



Focus™ Bright & True LED Light Technical Specs

Control Unit	Focus™ Xtend
Lamp weight	17 g
Lamp material	Aluminium
Operating time	6-24 hours depending on chosen light intensity
Charging time	Approx. 2.5 hours
LED diode	Focus™ Bright: Natural white (6,500 Kelvin) Focus™ True: True colour, low blue (5,000 Kelvin)
Maximum light intensity at 25 cm	Focus™ Bright: 42,000 lux Focus™ True: 40,000 lux
Light intensity at working distance 35 cm	Focus™ Bright: 19,000 lux (reduced glare) Focus™ True: 17,000 lux (reduced glare)
Light spot at working distance 35 cm	70 mm
Brightness setting	High, Medium and Low
Connection system	Focus™ Connect System
Cable type	Grey super flex
Warranty	3 years on Control Unit Xtend, Battery Xtend, Lamp unit True, Lamp unit Bright and & accessories 1 year on detachable cable and backplate