



Kepler Advanced Technical Specifications

Loupe lens system	Kepler, with higher light transmission than Kepler Kompakt
Lenses	Plano or Prescription
BlueLight Protection (BLP)	Yes (optional add-on)
Prisma Lens System	Yes (optional add-on)
Fine-focus available	Yes
Power-Up System	No, not needed (comes with 4 mags in one)
Magnifications	Turn'N'Click option of: 3.6x + 4.5x + 5.5x + 6.4x
Working distance	From 32 to 65 cm
Field of view (at 40 cm)	3.6x: 8 cm 4.5x: 7 cm 5.5x: 6 cm 6.4x: 5 cm
Depth of field (at age 50 years)	3.6x: 9.5 cm 4.5x: 9 cm 5.5x: 8.5 cm 6.4x: 5.5 cm
Mount system	TTL (Through the Lens)
Weight of two oculars	50.4 g
Loupe colour	Anthracite grey/Metallic grey
Loupe material	A Carbon enhanced composite housing and a titanium ring for changing magnification
Frames available	All