



# Kepler Reflekt

## Technical Specifications

Loupe lens system	Kepler Reflekt
Lenses	Plano or Prescription
BlueLight Protection (BLP)	Yes (optional add-on)
Prisma Lens System	Yes (optional add-on)
Fine-focus available	Yes
Magnification	4.0x
Working distance	From 35 to 57 cm
Field of View (at 45 cm)	7.2 cm
Depth of Field (at age 40 years)	9 cm
Mount system	TTL
Weight of two oculars with correction lens	36.9 g
Loupe colour	Black housing with titanium front
Loupe material	Magnesium and titanium
Frames available	All