



New Carbon v1 Be in your element.

Feather-light, yet strong and durable carbon-fibre frame.



Elevate your performance and style.

Superior comfort

Due to its ultra-low weight, it is easy to wear over long hours. The weight is only 50 g (Carbon v1 small with HD 2.3x loupes).

Durability

The perfect synergy of high-end materials. A carbon-fibre front and titanium temples for a unique combination of lightness and strength.

Unique look

The look of carbon-fibre is different from any other materials differentiating it not only in terms of technology but also in terms of design.

Anti-allergic

Carbon is a natural element (C is no. 6 in the Periodic system). It is non-poisonous and anti-allergic.




Frame specifications

Loupe range available	Kepler Advanced, Galilean HD & Essential
Frame material	Carbon-fibre
Temples material	Titanium
Weight	Only 50 g with HD2.3x (frame approx. 12 g)
Sizes	S, M & L
Suitable for Focus LED light system	Yes
Temples colours	Raw & Rose Gold



Are you now ready to experience true lightness?



Be in your element

Contact your local dealer to discover superior comfort, unrivalled lightness and a stylish design with your new Carbon v1 frame.

www.examvision.com

