



# PowerBand

## Technical Specs

Control Unit	PowerBand
Weight	195 g
Operating time	Varies with light intensity, 12-24 hours
Charging time	Approx. 2.5 hours
LED system	Can be used with lamp units Focus™ Bright and True and with Essential Nano
Brightness setting	High, Medium and Low
Connection system, output plug	MiniJack 3.5 mm
Input plug	Micro-B USB
Battery type	2x18650 lithium-ion rechargeable batteries with built-in overload protection
Battery voltage	3.6 VDC
Battery capacity	2x3450 mAh
Input	5 VDC/2A
Output (max)	2.9 VDC/400mA
Flicker-free	PowerBand is designed for ExamVision LED lights to provide flicker-free light.
Materials	Battery housing: PC-ABS V0 Neck strap: Silicone
Warranty	3 years