



EXAM **V**ISION

Extraordinary
performance
for medical use

Finest, custom-crafted
Danish loupes and
LED systems



Custom-crafted medical loupes for your needs



ExamVision offers versatile, all-day-wear tools to physicians, surgeons, and other medical professionals and students who require the ultimate optical loupe. Remarkably lightweight yet made from the most robust, high-end materials available.

Medical examinations and procedures require the use of superior optics and tools to facilitate freedom of movement, good posture and correct ergonomics throughout your working day.

Choosing the right loupe and magnification will optimise your performance in surgical and other medical procedures. Each solution is based on individual eyesight, need for magnification as well as physical working environment such as distance, tools and facilities available. So, your loupe is truly tailored and custom-crafted to you.

Your benefits in brief

All of our products are the result of a perfect synergy between technology and craftsmanship. We produce and assemble each loupe based on your individual measurements.

- **Excellent vision and precision**
Our optics help you achieve the crystal-clear image you need for high-quality, precision work.
- **Longer working life**
An improved ergonomic position via customised loupes is proven to avoid musculoskeletal problems.
- **Work satisfaction**
Achieved through better treatment results, faster work and happier, well-treated patients.
- **Style and functionality**
Five-time winner of the international product design award - Red Dot Design Award.
- **Peace of mind**
Five-year warranty on all of our loupes, with unrivalled dealer service support.
- **Customised solution**
Our specialists in optics and ergonomics ensure a solution tailored uniquely to you.

ExamVision loupe systems



reddot design award
winner 2010

Galilean HD

- Individual magnifications: 2.3, 2.8, 3.3, 3.8
- Perfectly balanced comfort
- Maximum field depth and width
- Ideal for all day wear
- Titanium oculars
- Any prescription possible



reddot design award
winner 2014

Kepler Kompakt

- Individual magnifications: 3.5, 4.6, 5.7
- Lightweight, ultra-compact design
- Suits the most demanding procedures
- Wide field of view, crystal-clear from edge to edge
- PowerUp available to upgrade magnification
- Lockable fine-focus ensures precision
- Freedom of movement via excellent depth of field



reddot design award
winner 2019

Kepler Advanced

- The only loupe with 4 magnifications: 3.6, 4.5, 5.5, 6.4
- Switch mag. easily via the Turn'N'Click system
- Very light and comfortable, for all day wear – 20 grams lighter than other comparable loupes
- Unparalleled sharpness from low to high magnification
- Made of titanium and carbon fibre-reinforced composite



ExamVision lights

Having the right light at the right spot is vital, both in terms of good vision and protecting your eyes from glare and harmful light. Loupes and light are interdependent and choosing the right headlight solution is almost as important, as choosing the optimum eye loupe solution. Your ExamVision dealer will guide you through personal consultation to ensure that you get the maximum benefit from your system.

Focus™ Bright

- Bright, neutral white $\leq 6,500$ Kelvin
- Light intensity 19,000 lux at 35 cm
- Light spot, 70 mm at 35 cm
- Three brightness settings
- Lamp weight 17 g
- Up to 6-24 hours operating time

Focus™ True

- Bright, neutral white $\leq 5,000$ Kelvin
- Light intensity 17,000 lux at 35 cm
- Light spot, 70 mm at 35 cm
- Three brightness settings
- Lamp weight 17 g
- Up to 6-24 hours operating time



Custom-crafted loupes and optics

ExamVision TTL loupes use the highest quality components and proactive coatings available, handcrafted to your exact specifications.

- Loupe angle up to 40° and custom working distance up to 60 cm, to improve posture and ergonomics.
- World class optics with the finest anti-reflective, protective coatings.
- Built-in prevention of chromatic aberration (colour separation).
- Carrier lens solutions – 100% customised: we offer Monofocal, Bifocal straight line and Multifocal/ Progressive lenses.
- Accommodates any prescription - as well as future prescriptions and changes to eyesight.
- Hi-Index lenses (lighter, thinner and more impact-resistant).
- BlueLight Protection (BLP), in both the ocular and carrier lenses available against harmful Blue-Violet light.
- Vision+Lens Protection (VLP) with built-in anti-reflection, filtering 99% of harmful UV light.

Frames

We strive to provide you with durable, yet fashionable frame options that will elevate your performance and complement your personal fashion preference. All ExamVision frames come with a 5-year warranty.

Panto v1



Sport v2



Carbon v1



Icon v2




Our Essence




At ExamVision we combine the rich heritage of Danish design with a commitment both to aesthetics and functionality. We come from Samsø, a tiny Danish island which is known worldwide for its renewable energy.

The peace and natural beauty of Samsø affords us the perfect environment in which to think big, work hard and pay meticulous attention to every single detail.

Try out ExamVision loupes for yourself and receive complimentary advice and guidance on vision, optics, protection and ergonomics at your ExamVision dealer.

Learn more about our magnifications



-  facebook.com/ExamVisionDK
-  instagram.com/examvision/
-  linkedin.com/company/examvision-aps

www.examvision.com

Get contacted

