

Kepler Kompakt and Kepler Advanced are direct view (prismatic) Kepler systems. These examples are to be used as a guide and represent a linear magnification, based on a working distance of 45 cm.



No magnification

Kepler Kompakt



Magnification	3.5x	4.6x	5.7x
Field of View	8 cm	7.5 cm	5.5 cm
Depth of Field up to	9 cm	6.5 cm	6 cm

Kepler Advanced



Magnification	3.6x	4.5x	5.5x	6.4x
Field of View	8 cm	7 cm	6 cm	5 cm
Depth of Field up to	9.5 cm	9 cm	8.5 cm	5.5 cm